***Instruction Manual
Manuale di istruzioni
Manuel d'instructions
Manual de instrucciones
Bedienungsanleitung***

ECO 16 Thermoreactor

F10100126, F10110126

**General Information / Informazioni Generali / Informations Générales / Información General /
Allgemeine Hinweise**

Before using the unit, please read the following instruction manual carefully.

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Do not dispose of this equipment as urban waste, in accordance with EEC directive 2002/96/CE.

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This unit must be used for laboratory applications only.

The manufacturer declines all responsibility for any use of the unit that does not comply with these instructions.

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This unit has been designed and manufactured in compliance with the following standards:

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Safety requirements for electrical equipment for measurement, control and for laboratory use

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IEC/EN 61010-1

Electrical equipment for laboratory use

UL 61010-1

General requirement - Canadian electrical code

CAN/CSA-C22.2 No.61010-1

VELP reserves the right to modify the characteristics of its products with the aim to constantly improving their quality.

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Contents / Indice / Index / Índice / Inhalt

1.	INTRODUCTION	5
2.	ASSEMBLY AND INSTALLATION.....	5
	2.1 ELECTRICAL CONNECTION.....	5
	2.2 START-UP	5
3.	OPERATING CONTROLS	5
	3.1 PROGRAMMING	5
	3.2 WORK CYCLE	5
	3.3 HEATING BLOCK TEMPERATURE VERIFICATION	5
	3.4 POWER INTERRUPTIONS	6
	3.5 SAFETY FEATURES	6
	3.6 END OF WORK CYCLE	6
4.	MAINTENANCE	6
	4.1 CLEANING	6
5.	TECHNICAL DATA	6
6.	ACCESSORIES / SPARE PARTS	6

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Thermo reaction can be carried out quickly and safely by placing the tubes into the thermoreactor heating block. Set the working temperature from room to 160 °C and the working time up to 199 minutes, or continuous. ECO 16 offers the possibility to thermoreact 14 tubes (Ø 16 mm) and 2 (Ø 22 mm) simultaneously.

2. Assembly and installation

Upon receipt and after having removed the packaging, please check the integrity of the instrument. The box includes:

- ECO 16
- Power supply cord
- Instruction manual

2.1 Electrical connection

Before connecting the instrument to the power supply, make sure that the values on the rating plate correspond to those of the power supply. Connect the unit to the power supply using the transformer supplied.

Ensure that the socket and the relative cut-off device conform to current safety norms and easy to reach.

2.2 Start-up

Use the on-off switch on the front of the unit to turn the unit on (position “-”). When switched on the unit all the leds blink for few seconds. When the leds are off the unit is ready to be programmed at different work temperatures and times.

3. Operating controls

Shift (A): it moves the messenger under all the parameters to be set during programming.

Increase (B): it shifts the value of the parameter corresponding to the messenger position.

Enter (C): it stores the values shown by the display.

Start (D): it sets going the work cycle.

Main switch (E): it turns on and off the instrument.

3.1 Programming

Turn on the main switch on control panel (E).

The unit performs in a few second a "Self-Test", after which: VELP SCIENTIFICA 20 °C *READY*

Pressing the Shift push-button the display shows two values contained in the memory (temperature and time):

Set	Temp.	Time
	100 °C	30'

Where:

Temp. (100 °C) is the temperature set. Range from room temperature to 160 °C.

Time (30') is the time of staying at the temperature set. Range from 0 to 199 minutes.

A continuous operation can be set by a 200, 300... value of time. A symbol === is displayed in this case.

The parameters set are stored by pressing Enter (C). When the unit is turned off the set values are retained in the memory.

3.2 Work cycle

It is set going by pressing the Start (D) and the display shows: Working 20 °C 60'

Where:

Working	Shows the operation of the work cycle.
20 °C	Is the present temperature of heating block.
60'	Is the time of staying at the set temperature.

When the set value of temperature is reached the value of time displayed starts to decrease with a resolution of one second.

At the end of work cycle a beeper warns and the display shows: *END*

When a new work cycle should start with a temperature of the heating block higher than the set value, the display shows:

VELP SCIENTIFICA 160 °C NOT READY

Where the 160 °C refers to the temperature of the heating block and NOT READY refers to the instrument operation.

If the Start push-button is pressed the cycle does not start and the display shows: NOT ACTIVE Temp. too high

Should a work cycle be interrupted the push buttons (C) and (B) must be pressed contemporaneously.

3.3 Heating block temperature verification

The temperature is controlled electronically which ensures stability without temperature oscillations around the set point.

The Pt 100 probe provides high precision and does not require calibration because the internal software performs a self-calibration every time the instrument is turned on. If, in spite of everything, a control is necessary, it is possible to measure

the temperature of the heating block with an external 3mm probe, introducing the 3mm probe pad into the hole in the upper part of the block

3.4 Power interruptions

The unit will register a black out or a lack of power supply as a manual stop. The work cycle will have to be restarted.

3.5 Safety features

The unit is protected against overheating by a thermostat that turns off the power supply when the temperature exceeds the safety temperature. Temperature setting is continuously controlled by electronic equipment that stops thermoregulation if anomalous measured are sensed. An alarm message appears on the display.

3.6 End of work cycle

At the end of the count-down the work cycle will stop automatically. To stop in every moment the work cycle turn off the general switch. When the work session is completed, should the unit remain connected to the power supply, switch it off using the on-off switch on the frontal panel.

4. Maintenance

No routine or extraordinary maintenance is necessary apart from periodically cleaning the unit as described in this manual. In compliance with the product guarantee law, repairs to our units must be carried out in our factory, unless previously agreed otherwise with local distributors. The instrument must be transported in its original packaging and any indications present on the original packaging must be followed (e.g. palletized).

It is the responsibility of the user, to properly decontaminate the unit in case of hazardous substances remaining on the surface or interior of the device. If in doubt about the compatibility of a cleaning or decontamination product, contact the manufacturer or distributor.

4.1 Cleaning

Disconnect the unit from the power supply and use a cloth dampened with an non-inflammable non-aggressive detergent.

5. Technical data

Power	550 W
Dimensions (WxHxD)	168x110x269 mm
Weight	3,8 Kg
Holes and diameter	14 (Ø16 mm) + 2 (Ø22 mm)
Construction material	Epoxy painted aluminum cast
Thermoregulation	Pid microprocessor
Selectable temperatures	room - 160 °C
Selectable times	0 - 199 min. - infinite
Set temperature reached	Visual signal
Count-down	Visual signal
End of cycle	Acoustic and visual signal
Overtemperature	On the block at 180 °C
Heating block temperature stability and homogeneity	± 0.5 °C
Temperature precision	± 1.0 °C
Enviromental ambient temperature	5 - 40 °C (41 - 104°F)
Storage temperature range	-10... + 60 °C (14 - 140°F)
Max humidity	85%
Pollution degree	2

6. Accessories / Spare parts

CA0091666	Test tube for sample decomposition 22 mm external diameter, NS 19/26 cone with glass cap
CA0091667	Condenser 200 mm type KS with 3 meters polyethylene tube
10000002	Absorption attachment for condenser
CM0091680	Set of 20 glass test tubes with internal Ø 14 mm and external Ø 16 mm
CA0091636	Holder for 12 round glass cells
A00001051	Safety shield
A00000362	Adaptor for probe Ø3mm for DK / DKL / ECO

10000667	ECO 16 Fuse 5x20 3A (230 V)
10000669	ECO 16 Fuse 5X20 5A (115 V)
10000232	Foot

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
We, the manufacturer VELP Scientifica, under our responsibility declare that the product is manufactured in conformity with the following standards:

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2011/65/UE (RoHS)

2012/19/UE (RAEE)

and satisfies the essential requirements of the following directives:

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- Machinery directive 2006/42/EC/ Macchine 2006/42/EC/ Machines 2006/42/EC/ Máquinas 2006/42/EC/ Maschinen 2006/42/EC

- Low voltage directive 2014/35/EU / Bassa tensione 2014/35/EU / Basse tension 2014/35/EU / Baja tensión 2014/35/EU / Niederspannung 2014/35/EU

- Electromagnetic compatibility directive 2014/30/EU/ Compatibilità elettromagnetica 2014/30/EU/ Compatibilité électromagnétique 2014/30/EU/ Compatibilidad electromagnética 2014/30/EU / Elektromagnetische Verträglichkeit 2014/30/EU

- plus modifications / più modifiche / plus modifications / más sucesivas modificaciones / in der jeweils gültigen Fassung.

Thank you for having chosen VELP!

Established in 1983, VELP is today one of the world's leading manufacturer of analytical instruments and laboratory equipment that has made an impact on the world-wide market with Italian products renowned for innovation, design and premium connectivity. VELP works according to **ISO 9001**, **ISO14001** and **OHSAS 18001** Quality System Certification.

Our instruments are manufactured in Italy according to the IEC 1010-1 and CE regulation.

Our product lines:

Analytical instruments

Elemental Analyzers
Digestion Units
Distillation Units
Solvent Extractors
Fiber Analyzers
Dietary Fiber Analyzers
Oxidation Stability Reactor
Consumables

Laboratory Equipment

Magnetic Stirrers
Heating Magnetic Stirrers
Heating Plates
Overhead stirrers
Vortex mixers
Homogenizers
COD Thermoreactors
BOD and Respirometers
Cooled Incubators
Flocculators
Overhead Shakers
Turbidimeter
Radiation Detector
Open Circulating Baths
Pumps

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VELP Scientifica Srl