***Instruction Manual  
Manuale di istruzioni  
Manuel d'instructions  
Manual de instrucciones  
Bedienungsanleitung***

**ECO 8 & ECO 25 Thermoreactors**

F101A0127, F101A0125

**General Information / Informazioni Generali / Informations Générales / Información General /  
Allgemeine Hinweise**

Before using the unit, please read the following instruction manual carefully.

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Do not dispose of this equipment as urban waste, in accordance with EEC directive 2002/96/CE.

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**This unit must be used for laboratory applications only.**

The manufacturer declines all responsibility for any use of the unit that does not comply with these instructions.

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**This unit has been designed and manufactured in compliance with the following standards:**

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Safety requirements for electrical equipment for measurement, control and for laboratory use

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**IEC/EN 61010-1**

Electrical equipment for laboratory use

**UL 61010-1**

General requirement - Canadian electrical code

**CAN/CSA-C22.2 No.61010-1**

VELP reserves the right to modify the characteristics of its products with the aim to constantly improving their quality.

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Thermo reaction can be carried out quickly and safely by placing the tubes into the thermoreactor heating block. The Start-Temp key offers the possibility to set 5 different work temperatures (70, 100, 120, 150 and 160 °C) whilst the Time key to set 4 different work times at the selected temperatures (30, 60, 120 or infinite). ECO 8 and ECO 25 offer the possibility to thermoreact respectively 8 and 25 tubes simultaneously (Ø 16 mm). ECO 8 model allows to hold also 1 test tube with Ø 22 mm.

## 2. Assembly and installation

Upon receipt and after having removed the packaging, please check the integrity of the instrument. The box includes:

- ECO 8 or ECO 25
- Power supply cord
- Safety cover
- Instruction manual

### 2.1 Electrical connection

Before connecting the instrument to the power supply, make sure that the values on the rating plate correspond to those of the power supply. Connect the unit to the power supply using the transformer supplied.

Ensure that the socket and the relative cut-off device conform to current safety norms and easy to reach.

### 2.2 Start-up

Use the on-off switch on the front of the unit to turn the unit on (position “-”). When switched on the unit all the leds blink for few seconds. When the leds are off the unit is ready to be programmed at different work temperatures and times.

## 3. Operating controls

|                        |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                 |
|------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>ON - OFF SWITCH</b> | It turns the unit on and off. If the switch is in the “0” position the unit is off; “-” means it is on.                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                         |
| <b>START-TEMP</b>      | When <b>Start-Temperature key</b> is pressed, a work cycle will start with the pre-set values and the corresponding leds will light: <b>Temperature 160 °C and Time 30 min.</b><br>The next pressing permits to select a different temperature: 160, 70, 100, 120 and 150 °C.<br>It's possible to select different work temperature until the heating block has not reached the set temperature and the count-down is not working. After that the key doesn't work anymore. If the selected temperature is lower than the running temperature of the heating block, (when the Start-Temp key is pressed) the work cycle will not start and all temperature leds will light intermittently, signalling an error. To correct the error, select a temperature that is higher than the running temperature or wait until the heating block cools down to a lower temperature than the temperature to be programmed. |
| <b>TIME</b>            | When the <b>Time key</b> is pressed it's possible to select the following different time: 120, 60, 30 minutes or ∞ (infinite). When a time is selected, the corresponding led will light. After switching on the instrument, the Time key permits to set the time at the set temperature but not to start the reaction cycle. The run time will start when the heating block reaches the set temperature.                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                       |
| <b>TEMP. LED</b>       | When a led lights it means that a work cycle is running at the selected temperature. If the leds are off a work cycle is not in progress.                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                       |
| <b>TIME LED</b>        | When a pre-set time is selected the corresponding led will light. When the led blinks it means that the selected temperature has been reached and the count-down has started.                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                   |

### 3.1 Heating block temperature verification

The temperature is controlled electronically which ensures stability without temperature oscillations around the set point. The Pt 100 probe provides high precision and does not require calibration because the internal software performs a self-calibration every time the instrument is turned on. Check the heating block temperature with a 3mm probe using the VELP Scientifica probe adaptor.

### 3.2 Power interruptions

The unit will register a black out or a lack of power supply as a manual stop. The work cycle will have to be restarted.

### 3.3 Safety features

The unit is protected against overheating by a thermostat that turns off the power supply when the temperature exceeds the safety temperature. Temperature setting is continuously controlled by electronic equipment that stops thermoregulation if anomalous measured are sensed. In case, all the leds start blinking and an acoustical warning signal sounds.

### 3.4 End of work cycle

At the end of the count-down the work cycle will stop automatically. To stop in every moment the work cycle turn off the general switch. When the work session is completed, should the unit remain connected to the power supply, switch it off using the on-off switch on the frontal panel.

## 4. Maintenance

No routine or extraordinary maintenance is necessary apart from periodically cleaning the unit as described in this manual. In compliance with the product guarantee law, repairs to our units must be carried out in our factory, unless previously agreed otherwise with local distributors. The instrument must be transported in its original packaging and any indications present on the original packaging must be followed (e.g. palletized).

It is the responsibility of the user, to properly decontaminate the unit in case of hazardous substances remaining on the surface or interior of the device. If in doubt about the compatibility of a cleaning or decontamination product, contact the manufacturer or distributor.

### 4.1 Cleaning

Disconnect the unit from the power supply and use a cloth dampened with an non-inflammable non-aggressive detergent.

## 5. Technical data

| Code Model                                          | F101A0127 ECO 8               | F101A0125 ECO 25         |
|-----------------------------------------------------|-------------------------------|--------------------------|
| Admitted power supply                               | 115 V - 230 V / 50–60 Hz      | 115 V - 230 V / 50–60 Hz |
| Power                                               | 135 W                         | 400 W                    |
| Dimensions (WxHxD)                                  | 135x95x230 mm                 | 155x95x275 mm            |
| Weight                                              | 2 Kg                          | 3,6 Kg                   |
| Holes and diameter                                  | 8 (Ø16 mm) + 1 (Ø22 mm)       | 25 (Ø16 mm)              |
| Construction material                               | Epoxy painted aluminum cast   |                          |
| Thermoregulation                                    | Pid microprocessor            |                          |
| Selectable temperatures                             | 70 - 100 - 120 - 150 - 160 °C |                          |
| Selectable times                                    | 30 - 60 - 120 min. - infinite |                          |
| Set temperature reached                             | Acoustic signal               |                          |
| Count-down                                          | Visual signal                 |                          |
| End of cycle                                        | Acoustic signal               |                          |
| Overtemperature                                     | On the block at 200 °C        |                          |
| Heating block temperature stability and homogeneity | ± 0.5 °C                      |                          |
| Temperature precision                               | ± 1.0 °C                      |                          |
| Time to reach 150 °C                                | 12 min.                       |                          |
| Enviromental ambient temperature                    | 5 - 40 °C (41 - 104°F)        |                          |
| Storage temperature range                           | -10... + 60 °C (14 - 140°F)   |                          |
| Max humidity                                        | 85%                           |                          |
| Pollution degree                                    | 2                             |                          |

## 6. Accessories / Spare parts

|           |                                                                                                  |
|-----------|--------------------------------------------------------------------------------------------------|
| CA0091666 | Test tube for sample decomposition 22 mm external diameter, NS 19/26 cone with glass cap (ECO 8) |
| CA0091667 | Condenser 200 mm type KS with 3 meters polyethylene tube (ECO 8)                                 |
| 10000002  | Absorption attachment for condenser (ECO 8)                                                      |
| A00001039 | Test tube extractor (ECO 25)                                                                     |
| CM0091680 | Set of 20 glass test tubes with internal Ø 14 mm and external Ø 16 mm (ECO 8 and ECO 25)         |
| CA0091636 | Holder for 12 round glass cells (ECO 8 and ECO 25)                                               |
| A00000362 | Adapter x probe Ø3mm for DK/DKL/ECO                                                              |

|          |                              |
|----------|------------------------------|
| 10000667 | ECO 8 Retarded fuse 5x20 3A  |
| 10000670 | ECO 25 Retarded fuse 5x20 8A |
| 10000232 | Foot                         |

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
We, the manufacturer VELP Scientifica, under our responsibility declare that the product is manufactured in conformity with the following standards:

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2011/65/EU (RoHS)

2002/96/CE (RAEE)

and satisfies the essential requirements of the following directives:

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- Machinery directive 2006/42/EC / Macchine 2006/42/CE / Machines 2006/42/CE / Máquinas 2006/42/CE / Maschinen 2006/42/EG

- Low voltage directive 2006/95/EC / Bassa tensione 2006/95/CE / Basse tension 2006/95/CE / Baja tensión 2006/95/CE / Niederspannung 2006/95/EG

- Electromagnetic compatibility directive 2004/108/EC / Compatibilità elettromagnetica 2004/108/CE / Compatibilité électromagnétique 2004/108/CE / Compatibilidad electromagnética 2004/108/CE / Elektromagnetische Verträglichkeit 2004/108/EG

- plus modifications / più modifiche / plus modifications / más sucesivas modificaciones / in der jeweils gültigen Fassung.

## Thank you for having chosen VELP!

Established in 1983, VELP is today one of the world's leading manufacturer of analytical instruments and laboratory equipment that has made an impact on the world-wide market with Italian products renowned for innovation, design and premium connectivity. VELP works according to **ISO 9001**, **ISO14001** and **OHSAS 18001** Quality System Certification.

Our instruments are manufactured in Italy according to the IEC 1010-1 and CE regulation.

Our product lines:

### Analytical instruments

Elemental Analyzers  
Digestion Units  
Distillation Units  
Solvent Extractors  
Fiber Analyzers  
Dietary Fiber Analyzers  
Oxidation Stability Reactor  
Consumables

### Laboratory Equipment

Magnetic Stirrers  
Heating Magnetic Stirrers  
Heating Plates  
Overhead stirrers  
Vortex mixers  
Homogenizers  
COD Thermoreactors  
BOD and Respirometers  
Cooled Incubators  
Flocculators  
Overhead Shakers  
Turbidimeter  
Radiation Detector  
Open Circulating Baths  
Pumps

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VELP Scientifica Srl