



# WIZARD

## Vortex Mixer

---

Exclusive and compact vortex mixer with settable speed, provided with the revolutionary IR sensor or continuous mode



# IR Sensor, New Conception of Vortex Mixers

Vortex mixer featuring the highly innovative and patented IR technology to reduce repetitive strains. IR sensor and Continuous mode available. Large choice of complementary accessories to accommodate tubes, plates and flasks.

- Vibration starts when tube is detected by IR (VELP patent)
- IR sensor and continuous modes
- Variable speed operation, electronically controlled up to 3000 rpm
- Optimum resistance to chemicals from technopolymer casing
- Excellent stability with zinc alloy base
- IP 42 protection against particles and liquids
- No fatigue as no pressure is required
- Eliminates repetitive strain injury



## INSTRUMENT - CODE

**WIZARD**

100-240 V / 50-60 Hz

**F202A0175**

\* Included: Mixing cup head

## TECHNICAL DATA

CONSTRUCTION MATERIAL	Zinc alloy and technopolymer
SUPPORT SYSTEM	3 anti-sliding feet
OPERATIONAL MODE	Sensor / Continuous
SPEED SETTING	Analog
SPEED CONTROL	Electronic
STIRRING SPEED	0 - 3000 rpm
ORBITAL DIAMETER	4,5 mm
TYPE OF MOVEMENT	Orbital
ELECTRONIC PROTECTION DEGREE CEI EN 60529	IP 42
POWER	15 W
WEIGHT	2.2 Kg (4.9 lb)
DIMENSIONS (WxHxD)	180x70x220 mm (7.1x2.8x8.7 in)



We reserve the right to make technical alterations  
We do not assume liability for errors in printing, typing or transmission

VELP Official Partner

## DESIGNED AND MANUFACTURED IN ITALY



**ITALY - HQ**  
Via Stazione 16  
20865 Usmate (MB) Italy  
Tel. +39 039 628811  
velpitalia@velp.com

**INDIA**  
velpindia@velp.com

**USA**  
155 Keyland Court, Bohemia  
NY 11716 - U.S.  
Tel. +1 631 573 6002  
velpusa@velp.com

**CHINA**  
Xinlong Rd Building 28, Lane 1333  
Shanghai city - China  
Tel. +86 180 1755 7329  
velpchina@velp.com