

VELP Scientifica Glass Crucibles



Glass crucible P0
6pcs/box

Glass crucible P1
6pcs/box

Glass crucible P2
6pcs/box

Glass crucibles P3
6pcs/box

FILTER	P0	P1	P2	P3
NOMINAL PORE	Size 160-250 micron	Size 100-160 micron	Size 40-100 micron	Size 16-40 micron
EQUIPMENT	FIWE – CSF6	FIWE – CSF6	FIWE – CSF6	FIWE – CSF6
EXAMPLE OF APPLICATIONS FIELDS	Special crucible; Filtration of very coarse precipitates	Special crucible; Filtration of coarse precipitates	Standard crucible; preparatory work with crystalline precipitates	Special crucible; analytical work with medium-fine precipitates
CODE	A00000138	A00000139	A00000140	A00000137



www.velp.com