

VELP Scientifica Automatic Digestion Units



DKL 8



DKL 12



DKL 20



DKL 42/26

| | | | | | |
|--|-------------------------|------------------------------------|------------------------------------|------------------------------------|------------------------------------|
| Number of positions | | 8 x 250 ml | 12 x 250/400 ml | 20 x 250 ml | 42 x 100 ml |
| Lift movement | | Automatic | Automatic | Automatic | Automatic |
| Temperature range | | Ambient to 450 °C | Ambient to 450 °C | Ambient to 450 °C | Ambient to 450 °C |
| Stability and precision of heating block temperature | | ± 0.5 °C | ± 0.5 °C | ± 0.5 °C | ± 0.5 °C |
| Program library | | 30 standard + 24 user-programmable | 30 standard + 24 user-programmable | 30 standard + 24 user-programmable | 30 standard + 24 user-programmable |
| Selectable ramps | | up to 4 per program | up to 4 per program | up to 4 per program | up to 4 per program |
| Digestion time range | | from 1 to 999 minutes | from 1 to 999 minutes | from 1 to 999 minutes | from 1 to 999 minutes |
| Supplied with | Lift | • | • | • | • |
| | Suction cap & drip tray | • | • | • | • |
| | Sample rack | • | • | • | • |
| | Test tubes | • | • | • | • |
| Dimensions (WxHxD) mm | | 210x690x540 | 266x690x540 | 322x690x584 | 322x690x584 |
| | (in) | (8.3x27.2x21.3) | (10.5x27.2x21.3) | (12.7x27.2x23.0) | (12.7x27.2x23.0) |
| Weight | Kg | 19.7 | 23.3 | 30.8 | 33.5 |
| | (lb) | (43.5) | (51.4) | (68.0) | (74.0) |
| Power Supply | | 115 or 230 V | 115 or 230 V | 230 V | 230 V |
| Power | | 1150 W | 1500 W | 2300 W | 2300 W |



www.velp.com